Carsuh Acuproductos y Afines, S.L.

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Suhpresor®, the active foot reflexology device from **Suh Seung Yeul**.

Thanks to its exclusive and patented design, with its more than 300 pressure points engraved in it, you will be able to initiate your own personalized foot reflexology treatment in a few steps. Read these instructions carefully, and consider this information before using your Suhpresor® for the first time in order to obtain optimal results and avoid injury. If you have any doubts about it, please get in touch with the manufacturer and/or distributor.

## WHAT IS REFLEXOLOGY?

Reflexology is a technique based on the same principles as acupuncture. It consists of the application of pressure over specific areas of the body for one's health benefit. The most popular branch is foot reflexology. Discomfort, pain, even illness may occur when an imbalance is produced inside one of our organs or systems. Through its nerve endings, different parts of our organism have their reflection in our feet [See figure 4]. Using foot reflexology, by applying pressure in specific areas, one can stimulate the autonomous nervous system, thus creating a positive impact in our circulatory, lymphatic and nervous systems, among others.

## WHAT IS SUHPRESOR AND WHAT IS IT FOR?

Suhpresor® is the final result from years of research and development conducted by the acupuncture and acupressure professional **Suh Seung Yeul**, with experience since 1972.

When used correctly, one can apply an active foot reflexology pressure massage that promotes physical and mental health, helping with different ailments such as anxiety, stress, back pain, headache, digestive or menstrual disorders, insomnia, hormonal imbalance, etc. The majority of existing devices in today's market are designed for a passive use (seated, laying down or applied by others), instead of a personalized technique. Suhpresor®, with its exclusive and unique design, makes reflexology become an active and unique experience: it is you who exerts pressure over specific areas of your foot with your body weight and movement. This innovative device has been designed to adapt in the most efficient way to modern life's needs. It is for people who do not have time to go see a reflexologist due to work or personal reasons. It is compact, light and easy to transport. It's always ready to use whenever you so desire.

A great way to help with different diseases, disorders and energetic imbalances of the organism is to have a correct, active and periodic activity with Suhpresor®, accompanied by the exercises also designed by the renowned acupuncturist **Suh Seung Yeul**.

## HOW TO USE SUHPRESOR

Regardless of the specific use you may wish to give to Suhpresor to treat a particular pain or individual problem (see figures 5, 6, 7 and 8), it is fundamental to always apply the **basic essential massage** beforehand. Its duration must be twice the time allotted for the specific

treatment you want to apply.

### ESSENTIAL BASIC MASSAGE:

It is the basic and general exercise for all ailments that must be applied before each specific treatment. It consists of a general massage using Suhpresor® across the entire foot sole on both feet, following this example (see figures 1, 2 and 3).

Both for the essential basic massage as well as the specific massage, place yourself in the position that best adapts to your needs and ability (see example from figure 1a). This way, supporting yourself with your hands, you will be able to compensate the pressure your body applies to your feet at your best convenience.

Once your feet are adapted to the pressure massage, the most recommended position is standing in upright position (see figure 1B).

Place the sole of your feet on top of the device. Start stepping on it, moving your feet with short steps forward and backwards (from the tips of your toes to your heels), and from left to right (see figure 2): this way you are applying a pressure massage to the entirety of your soles and lateral points on both feet (saphenous and sural). The corrugated waves of the device at the center and on the sides work towards this objective (see figure 3).

### RECOMMENDED TIME OF DAILY USE:

As a general rule, and only when your feet are used to it, the recommended time for daily use is a minimum of 10 minutes, and a maximum of 30 minutes per day and person. You may wish to spread out that time along different parts of the day and in different sessions. To obtain the best results with Suhpresor®, it is recommended to use it in two daily sessions. It is also advisable to use Suhpresor® with bare feet, or at least, the thinnest piece of clothing possible.

## CORRECT USE OF SUHPRESOR

At first, and until your feet are familiarized with the pressure massage, it is possible that some people may experience some discomfort, especially those overweight or with more sensibility in their feet. It is recommended during first use to start placing your foot gently, in short sessions, 2 or 3 minutes long, and rest. Afterwards you may repeat as many times as necessary (without exceeding the 30 minutes of daily recommended use). Once your feet are familiar with the massage, gradually increase the daily use of Suhpresor® until you reach the amount of time that fits your needs.

Please remember: The use of Suhpresor® must not exceed the recommended amount of 30 minutes per day per person.

Use the reflexology exercises only on your feet soles, with Suhpresor® always placed at foot level, over a flat, stable, non-slipping surface. In order to keep a vertical position and not lose balance and avoid accidents or injury caused by possible falls, grab hold of a wall, piece of furniture or any other element in your surroundings that is firmly fixed on the floor (see figure 1a).

## MAINTENANCE

Each Suhpresor® has been manufactured using solid hardwood treated with a natural water varnish, to provide maximum quality, hygiene and product durability. Nonetheless, with use, it is expected that dirt will accumulate on its

surface. To clean it we only recommend the use of a soft brush, soap and water. Rinse afterwards with clean water, dry it with a cloth and let it dry in an open and dry area.

## CONTRAINDICATIONS AND WARNINGS

The use of Suhpresor® does not have known secondary effects, but it is not recommended for use in these cases: Presence of wounds, feet ulcers, clotting disorders, medical history of pathological fractures, diabetic foot, people with fractures or sprained ankles that have not yet fully healed, and people with reduced mobility. Suhpresor® must not be used by underage people without the appropriate adult supervision. If you have any doubt, check with your physiotherapist, reflexologist or doctor. To guarantee the hygienic conditions of the product, it is necessary to follow the cleaning indications described in the MAINTENANCE section, for each new user or/and after a prolonged amount of time without use. Suhpresor® must always be used in flat, level, non-slippery surfaces, free of obstacles around it.

Suhpresor® must be kept in a dry area, far from contamination sources of any kind.

It is forbidden to use deteriorated or broken devices that may present any type of risk to its users.

Suhpresor® has been designed and manufactured only and exclusively for the uses indicated in this brochure. It is forbidden to use for any other activity not described here.

For more information, consult with your distributor and/or manufacturer, and visit [www.suhpresor.com](http://www.suhpresor.com)

Once Suhpresor® has finished its life cycle, it must be deposited in the Urban Residues Container.

As far as the brochure and box, they must be deposited in the blue container (Paper, cardboard).

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**Insomnia, stress and anxiety**  
 Please abide by the basic and general usage instructions and dedicate a quarter (1/4) of the total usage time to step on the following highlighted area.

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**Digestive, constipation and diarrhea**  
 Please abide by the basic and general usage instructions and dedicate a quarter (1/4) of the total usage time to step on the following highlighted area:

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**Migraines and headaches**  
 Please abide by the basic and general usage instructions and dedicate a quarter (1/4) of the total usage time to step on the following highlighted area:

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**Rheumatism and blood circulation**

Please abide by the basic and general usage instructions and dedicate a quarter (1/4) of the total usage time to step on the following highlighted area.

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## HELPS, AMONG OTHERS, IN:

- Nervous system
- Anxiety and depression
- Sleep and rest
- Leg and foot inflammation
- Hypertension
- Cholesterol
- Constipation
- Menstrual cramps
- Headache and Cephalea
- Improvement of blood flow
- Fluids circulation
- Obesity
- Benefits for all organs and viscera

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